#### Advances in quantitative Rietveld Analysis XRPD for Minerals and Mining Applications



Dipl. Min. Alexander Seyfarth BRUKER AXS Inc., Madison WI, USA

Dr. Arnt Kern BRUKER AXS, Karlsruhe, Germany



#### Outline



- Why X-Ray Powder Diffraction (XRPD)?
- XRPD application areas and capabilities
- Recent advances in quantitative phase analysis with XRPD
- Example applications
  - Process / production control
- Conclusions



- X-Ray Powder Diffraction (XRPD) is an analytical tool for materials characterization, including but not limited to
  - qualitative phase analysis (phase identification),
  - quantitative phase analysis,
  - crystal structure determination and refinement
  - and much more
- XRPD is sensitive to the <u>crystal structure</u> of each phase present in the sample
- The emphasis of this presentation is on <u>quantitative phase</u> <u>analysis</u>



Alternative (non XRPD) methods for quantitative phase analysis:

- <u>Point counting</u> using an optical microscope, scanning electron microscope or electron microprobe, now usually combined with digital image analysis
  - rather slow, difficult on-line automation->not usable for process
  - surface sensitive, can result in poor statistics
  - limited by fine grain size

# Modern "automated" or "quantitative" mineralogy



Image processing (EDS and BSE) to obtain:



- Chemical assay(s)
- Modal mineral proportions

• Grain size

+

- Mineral associations and liberations
- Porosity

 $\rightarrow$  important information for mining and processing

#### Ultra fast element mapping + Mineral ID Automated MINERALOGY (SEM/BSE)





- Mineral
- 20 kV / 10 nA
- 250 kcps
- 1024 x 768
- 15 min (1 detector)





Example: TiO<sub>2</sub>



Anatase (with calcite)



**Rutile (with hematite)** 

"Give rutile and anatase to chemists and they will tell you they are both 100% TiO<sub>2</sub>" (lan Madsen, CSIRO)







Example: Iron, iron oxides, iron hydroxides



### X-Ray Diffraction and Scattering



= {

ΞI



- Peak positions and intensities are functions of the crystal structure of a crystalline phase
- In mixtures, intensities are related to phase abundance
  - ⇒ Quantitative phase analysis
- An X-ray powder pattern is characteristic for a crystalline phase with its particular elemental composition <u>and</u> crystal structure
  - ⇒ "Fingerprint" phase identification



Why is this important?





- Materials properties are not solely determined by their chemistry as e.g. determined by XRF, but by its mineralogy, i.e. the crystal structure (s) of the constituent compound (s)
- Crystal structure governs properties such as
  - Crystal habit / morphology
    - Crystal surfaces
  - Surface charge distribution
    - Hardness
    - Density
      - ...

- Grindability
- Flowability
- Solubility
- Floatation properties

• ...

Knowledge of these properties is a prerequisite for optimum processing

### XRPD application areas and capabilities

#### **XRPD** capabilities

- Sample amounts: μ–grams (micro-diffraction) up to grams
- Ideal grain size required: <10μ
  - XRPD is sensitive to fine grain size
  - XRPD cannot provide information about particle size and shape, mineral association and liberation
- Quantitative analysis:

4/2/2011

- linear concentration range from 0.1-3%\*) to 100%
- typical accuracy 0.1-3%\*) and reproducibility <0.1%\*) absolute
- typical detection limits: 0.1-1%\*)
- \*) depends strongly on sample presentation and sample properties, such as elemental composition (scattering power), crystal structure symmetry, degree of crystallinity









Recent advances in quantitative phase analysis with XRPD



#### The Rietveld method

- The Rietveld method generates a calculated diffraction pattern that is compared with the observed data
  - Qualitative phase analysis required
- The differences between observed and calculated diffraction patterns are minimized using least-squares procedures



# Quantitative X-ray Mineralogy using the Rietveld Method

- Hill & Howardt (1987)
  J. Appl. Cryst. 20, 467-74
  - all phases identified
  - all phases crystalline
  - all crystal structures known
- Benefits
  - No need for artificial calibration mixtures
  - ZMV is the calibration constant
  - ZMV known from crystal structure

 $\frac{S_{\alpha}\left(ZMV\right)_{\alpha}}{\sum_{k}^{n}S_{k}\left(ZMV\right)_{k}}$  $W_{\alpha} =$ 

- W= Weight %
- S = Rietveld scale factor
- Z = No. of formula units in unit cell
- M = Molecular mass of formula unit
- V = Unit cell volume

Recent advances in quantitative phase analysis with XRPD



A new generation of Rietveld software: TOPAS (since 1997)

In addition to RECIPE based ONE button use:

- A convolution based <u>instrument function approach</u> for describing observed X-ray line profile shapes
   Fundamental Parameters Approach
  - REDUCES PARAMETERS which need to be fitted DRAMATICALLY and enables REAL STABLE REPEATABLE refinements
  - No divergence approach....

Recent advances in quantitative phase analysis with XRPD



Fundamental Parameters Approach

- The observed line profile shapes in a powder pattern are calculated from the <u>known instrument geometry</u>
- This allows a more reliable decomposition of peak overlaps at much higher degrees of peak overlap, compared to traditional analytical profile functions (e.g. pseudo-Voigt, PearsonVII)
- The number of refineable profile parameters and therefore parameter correlation is significantly reduced
  - Analytical profile fitting : ~7 (U,V,W,X,Y,Z,Asymmetry)
  - Fundamental parameters approach : ~1-2 (size, strain)
- Complex line profile shapes as found in clays can be modeled

#### Quantitative X-ray Mineralogy Amorphous material – PONKCS

- Use if spiking not feasible
- Quantification of Phases Of No Known Crystal Structure
  - amorphous
  - unknown or partly known structure
- Scarlett & Madsen (2006)Powder
  Diffraction 21(4), 278 284
- Calibration of an unknown phase α via internal standard s in TOPAS as
  - unindexed peaks phase
  - indexed hkl-phase

- S = Rietveld scale factor
- $Z \neq$  No. of formula units in unit cell

 $\left(ZMV\right)_{\alpha} = \frac{W_{\alpha}}{W_{\alpha}} \cdot \frac{S_{s}}{S} \cdot \left(ZMV\right)_{s}$ 

- M = Molecular mass of formula unit
- $V \neq$  Unit cell volume
  - ZMV known for standard, but calibrated for unknown





#### Example applications: production control / industrial





#### Quantitative Rietveld Analysis Typical Recent Subjects



- Mineral processing products at mining operations
- Composition of ores
- Mine tailings and waste rocks (environmental mineralogy)
  - Acid mine and rock drainage
- Mineralogy of asbestos mine tailings (CO<sub>2</sub> sequestration)
- Undesirable deposits, clogs, etc. in furnaces, boilers, pipes, etc.
  - E.g. early warning of blockage
- Mineralogy of (exotic) slags
- Miscellaneous corrosion products

#### Quantitative Rietveld Analysis Porphyry Copper Deposit Host Rock



• E.g. determination of the acid producing / neutralisation potential



#### Quantitative Rietveld Analysis Gold Mine Waste Rock



• E.g. determination of the acid producing / neutralisation potential



# Objective: XRD to support mine operation and ore processing







- Provide data for better planning and forecasting
  - resources
  - estimate ore reserves
  - ore control
  - haulage
  - energy / chemicals consumption

#### Quantitative Rietveld Analysis Phosphate Ore





#### Quantitative Rietveld Analysis Phosphate Ore



Phase	As extracted	Concentrate	Tailing
• Quartz	33.10	2.05	41.81
• Hematite	1.42		5.54
• Hydroxyapatite	38.20	89.20	17.55
• Dolomite	3.33	2.40	2.59
• Calcite	2.93	4.40	2.40
Goethite	6.23		9.86
Vermiculite	4.83		8.14
• Ilmenite	3.49		1.48
• Anatase	1.31		1.70
• Barite	0.41	1.95	0.64
• Diopside	2.36		3.75
Microcline	2.40		4.43

#### Quantitative X-ray Mineralogy Results reconciliation



- Accuracy of the XRD method can be validated by comparing against independent methods
  - chemical analysis (XRF, ICP-MS, AA, ...)
  - optical microscopy (quantitative point counting)
  - SEM-EDS

XRD and QUEMSCAN secondary minerals list



Quantitative Rietveld Analysis CO<sub>2</sub> Sequestration



- Wilson, S.A., Raudsepp, M. & Dipple, G.M. (2006).
  American Mineralogist 91, 1331-1341.
- The sequestration of anthropogenic CO2 may be required to meet Canada's commitment to the Kyoto Protocol
- The <u>carbonation of serpentine-group minerals</u> in ultramafic mine tailings presents an opportunity to implement carbon sequestration in the mining industry
  - Globally, ultramafic mines could sequester 10<sup>8</sup> tonnes of CO<sub>2</sub>/year
- The trouble with serpentine: Stacking disorder no reliable crystal structure

#### Quantitative Rietveld Analysis Structureless Fitting of Crysotile





#### CO<sub>2</sub> Sequestration Clinton Creek Mine, BC, Canada









- XRPD is a direct and accurate analytical method for determining the presence and absolute amounts of mineral species in a sample
- STANDARDLESS QUANTIFICATION is a reality and PROCESS ready with accuracies of better than 5% relative

 Significant advances have been made in quantification of disordered materials (e.g. clays) with new functionality in TOPAS

• IN LEACHING (HEAP) XRD is used as THE CONTROL TOOL

#### Conclusion



#### • Important limitations are

- The relatively high lower limit of detection, particularly for poorly crystalline phases
- The requirement for <u>appropriate sample preparation</u>. Overgrinding may destroy soft phases, resulting in an underestimation.

Note, that microscopes/microprobes, XRF, and XRD are highly complementary methods, but they look at different samples!





## **QUESTIONS**?

#### References



Scarlett, N.V.Y. & Madsen, I.C. (2006) *Quantification of phases with partial or no known crystal structure* Powder Diffraction, 21(4), 278-284

Wilson, S.A., Raudsepp, M. & Dipple, G.M. (2006). American Mineralogist 91, 1331-1341.

Webinars (recorded) or Online Classes

TOPAS WEBINAR, 2008 , please contact Karen.Roscoe@bruker-axs.com

CEMENT WEBINAR, 2008, please contact Karen.Roscoe@bruker-axs.com



www.bruker.com

© Copyright Bruker Corporation. All rights reserved